


‘BREAK GLASS’ FUNCTION

IntelConnect and IntelViewer

This function allows access to patient images when you are not the primary referrer.

1. Select **Break-Glass Access**

IntelConnect:




IntelViewer:




2. Acknowledge and accept terms of the Confidentiality Agreement

3. Search for the patient using **surname, first name** and **DOB (YYYYMMDD)**



4. Click **Access Patient Details**
5. You will then be returned to the 'Search' screen. Select **Search** to view all studies for the patient



To open patient file click anywhere in the patient line

To open images click anywhere on the image